

China Dental Lab Acrylic Resin Temporary Crown Dentures Material

Our Product Introduction

for more products please visit us on adschinadentallab.com

Basic Information



Product Specification

- Turnaround Time: 3 Days In Lab
- Warranty: 5 Years
- Shipping Way: UPS (preferred), Fedex And DHL Linex. Armex Are Also Available
- Payment: T/T, HSBC, PayPal, MoneyGram, Western Union
- Highlight: China Dental Lab Acrylic Resin Temporary Crown , Acrylic Resin Temporary Crown Material, China Dental Lab acrylic resin dentures



More Images



Product Description

Description

Acrylic Resin Temporary Crown: The Ultimate Solution for Temporary Dental Crowns

Are you looking for a reliable and effective solution for temporary dental crowns? Look no further than our Acrylic Resin Temporary Crown! As a leading manufacturer in China, we understand the importance of delivering high-quality dental products to our clients worldwide.

Our Acrylic Resin Temporary Crown is a superior solution that offers durability, strength and versatility to dental practitioners. Allow us to introduce you to the benefits of our product.

1. Excellent Aesthetic Appearance

Our Acrylic Resin Temporary Crown is designed to match the natural appearance of teeth. With a translucency and shine that mimics the appearance of natural enamel, this product provides an enhanced cosmetic look that patients are sure to approve of.

2. Durability and Strength

Our Acrylic Resin Temporary Crown is engineered with high-quality materials to provide maximum strength and durability. It is resistant to fracture, chipping and other damage, ensuring long-lasting functionality for patients.

3. Easy to Use

Our Acrylic Resin Temporary Crown is designed to make dental procedures quick, easy and efficient. It is easily shaped to the contours of the teeth and can be trimmed and polished to achieve the desired size, shape and feel.

4. Affordable and Cost-Effective

Choosing our Acrylic Resin Temporary Crown means you can benefit from an affordable solution without compromising on quality. It represents an excellent value for cost-conscious dental practitioners who do not want to compromise on quality.

5. Versatility

We understand that dental procedures can vary greatly, which is why our Acrylic Resin Temporary Crown is highly versatile. It can be used for a wide range of procedures, including to protect the tooth after a root canal, to hold a dental bridge in place, or as a temporary measure for dental implants.

6. Safety and Comfort

Our Acrylic Resin Temporary Crown is designed with patient comfort and safety in mind. It is minimally invasive, lightweight and comfortable to wear while providing necessary protection for the teeth.

In summary, our Acrylic Resin Temporary Crown is a comprehensive solution that offers excellent aesthetics, durability, and versatility at an affordable price point. It is easy to use, comfortable, and highly effective at protecting teeth during a variety of dental procedures.

As a trusted manufacturer based in China, we are committed to delivering the highest quality products to our clients worldwide. We use only the finest materials and advanced manufacturing techniques to produce dental products that exceed industry standards.

We invite you to partner with us and experience the advantages of our Acrylic Resin Temporary Crown today!

Production Introduction:

Polymethyl methacrylate acrylic (PMMA) milled crowns are an adequate, temporary, fast and reliable way to save chairside time, while the hardness and resistance of the material is perfect for the preparation of provisional pieces, designed for use until a definitive restoration can be made.

China Digital Dental lab fabricated provisional restorations are milled from PMMA (Polymethylmetacrylate). Our PMMA provisionals are CAD/CAM generated and can be created from either conventional impression or a digital IOS workflow.

Product Parameter

Customized

Product Feature and Application

These PMMA temporary crowns or bridges are recommended for long-term provisional needs such as implant or surgical cases. They are more durable than conventional provisionals and can be fabricated to receive titanium cylinders for use in the provisionalization of screw-retained implant cases.

Key Characteristics

- Polymethyl methacrylate material
- Biocompatible
- Highly resistant to abrasion
- CAD/CAM manufactured – for accuracy of fit
- A long-term temporary option, reducing chair side time

Deliver,Shipping and Serving

We accept STL files.

We accept normal impressions, you can send us by Fedex,UPS,DHL or Linex

We finish the jobs in 4-5 days once we receive your cases

We accept TT.Paypal





Why choose us?

Our Metal Free Crowns products are manufactured with a focus on safety and reliability.

We repay our customers' trust with high quality enterprise management, excellent Acrylic Resin Temporary Crown and perfect after-sales service.

Our Metal Free Crowns products provide patients with a long-lasting and reliable solution.

We believe that excellent service quality and good customer relationship are the basis for the survival and development of our company.

Our Metal Free Crowns products are designed to provide patients with a long-lasting solution without sacrificing aesthetics.

Over the years, professional companies from all over the world take us as their long-term and stable partners.

Our company has a proven track record of supplying Metal Free Crowns products to dental professionals worldwide.

Over the years, our company has been upholding the fundamental of excellent quality, technological progress as the core, green environmental protection as the responsibility, beyond innovation as the goal.

We place a strong emphasis on product innovation to offer patients the best Metal Free Crowns solutions.

Our company relies on its own professional technology, quality casting enterprise integrity, perfect and thoughtful after-sales service occupies the trust of the majority of customers.

Hot Tags: acrylic resin temporary crown, China acrylic resin temporary crown manufacturers, suppliers, factory, zirconia full arch bridge, palatal expander loose teeth, orthodontic mouth widener, palate expander on bottom, orthodontic spacer appliance, bonded expander orthodontics

ADS DENTAL LAB 卓艺美学中心 **ADS Dental Lab**

+8613430498390

info@chinaadsdentallab.com

adschinadentallab.com

Building 4, Hengguangyao Industrial Park, Yonghe Road, Fuyong, Shenzhen, China 518103